



July/August 2021

**3<sup>rd</sup> Ocular Camera & Telepathology Software  
Specifications, Pricing, Terms & Conditions**

Special  
**Special**  
Proposal

Send Email with PO to: [info@augmentiqs.com](mailto:info@augmentiqs.com)

## Camera & Software Specifications:

Resolution	20 MP or 5 MP (4K+ equivalent)
Sensor size	1 " diagonal
FPS	Up to 75
Shutter	20mp Rolling Shutter / 5mp Global Shutter
Bit depth	Up to 12 bit
Interface	USB 3.0
Software	<p>License includes software updates.</p> <p>Image capture (TIFF, PNG, JPG, BMP), Annotations, Morphometrics, Database, Telepathology, QuPath algorithms Ki-67 &amp; WCD</p> <p>Note: 5MP camera version supports manual scanning software.</p>

## Price List:

Line Item	Euro
3 <sup>rd</sup> Ocular Camera	3,300
Annual Telepathology License	2,400
SW Setup Support	400*
Morphometric Calibration Slide	200
Shipping	200
<b>Total Price</b>	<b>6,500</b>
<b>Discount</b>	<b>-3,000</b>
<b>Discounted Price</b>	<b>3,500</b>

\* Augmentiqs will assist in setting up the telepathology software after the user or microscope service provider installs the camera.

## Optional Accessories:

Line Item	Euro
Calibration Certificate for Morphometric Slide	150
Color calibration kit including color calibration slide	500
SOP documentation	750
Sending images to SMS, one year license	150
Mindpeak AI Algorithms Ki-67, ER/PR, HER2, one year license	Ask for pricing
Integration of other 3 <sup>rd</sup> party software/AI	Ask for pricing
Aira Matrix Algorithms	Ask for pricing
Manual WSI application	Ask for pricing

## Relevant Links for Additional Information:

Augmentiqs	<a href="http://www.augmentiqs.com">www.augmentiqs.com</a>
Mindpeak	<a href="http://www.mindpeak.ai">www.mindpeak.ai</a>
Aira Matrix	<a href="http://www.airamatrix.com">www.airamatrix.com</a>
Microvisioneer	<a href="http://www.microvisioneer.com">www.microvisioneer.com</a>

## SALES AGREEMENT TERMS & CONDITIONS

Augmentiqs Medical Ltd. (“Company”)

### ***General***

- All prices are in Euro, excluding VAT. Payment can be made in USD based on the current exchange rate
- The agreement does not include any purchase or license to use any of the Company’s proprietary rights or intellectual property unless specifically determined otherwise.
- This price quote is valid until August 31, 2021, or until end of supply. Company reserves the right to terminate this offer at any time.
- Company will ship Products to Buyer within 14 days of Payment.

### ***Limited Warranty***

- Company warrants that all products will be free from defects in material and workmanship, under normal use and service, for a period of six (6) months after the date of shipment.
- Warranty shall become void in the following cases: (i) misuse of the Products or any part thereof not in accordance with the Company’s written or verbal instructions, or (ii) any mishandling, tampering, opening or an unauthorized attempt to fix, or modify the Product or any part thereof, or (iii) accidents, act of war or terror, force majeure, etc.

### ***Return of Merchandise***

- Buyer may return merchandise for a partial refund within 30 days of receipt of goods. Refund includes only the cost of goods, but does not include shipping. 10% restocking fee of hardware costs will be deducted, and cost of installation support.
- Policy will not apply in event of (i) misuse of the Products or any part thereof not in accordance with the Company’s written or verbal instructions, or (ii) any mishandling, tampering, opening or an unauthorized attempt to fix, modify or installation of the Product or any part thereof, (iii) changes in Buyer IT policy that prevent use of the Product and/or one of its software modules, (iv) failure to provide

PC and/or Internet firewall settings that meet the recommended requirements listed in this PO agreement or (iv) accidents, act of war or terror, force majeure, etc.

### ***Liability***

- UNDER NO CIRCUMSTANCES WILL THE COMPANY BE LIABLE FOR DIRECT, INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES RELATING TO USE OF THE PRODUCT, INCLUDING, WITHOUT LIMITATION, LOSS OF PROFIT, LOSS OF USE, LOSS OF REVENUE, OR DAMAGES TO PERSON, BUSINESS OR REPUTATION ARISING FROM USING THE PRODUCTS, OR THE BREACH BY BUYER OF ANY TERM, REPRESENTATION OR WARRANTY OF THIS PURCHASE ORDER, OR ANY OTHER OBLIGATION TOWARD ANY THIRD PARTY AT LAW OR IN EQUITY, WHETHER OR NOT COMPANY WILL HAVE BEEN MADE AWARE OF SUCH LOSS OR DAMAGE.
- Augmentiqs is a microscope accessory, and is not intended for diagnostic use.

### ***User Hardware Requirements***

- The Buyer is solely responsible for assuring microscope and PC meet the recommended minimum operating requirements.
- Microscope integration – Buyer is solely responsible for installing the camera on their microscope.
- ***PC: Augmentiqs software is installed on BUYER'S existing PC, Windows 7 or higher, 64 bit (hi-capacity SSD storage with high throughput).***
  - CPU: i7-6770HQ or more capable (4 cores, 8 threads)
  - RAM: 8 GB DDR3 (for AI applications we recommend 16GB)
  - HDD: 100GB
  - Display: Full HD (1920 by 1080-pixel resolution) display with IPS technology for proper color gamut support
  - Ports: USB 3.0 (X1),
  - Internet upload speed of at least 2MB/s.
  - Internet-facing firewall should allow outbound and inbound TCP ports 443.